

CHILL TOWEL

DESCRIPTION

THORZT COOLING APPAREL IS DESIGNED TO ASSIST WITH HEAT STRESS PREVENTION BY PROMOTING THE BODY'S THERMOREGULATION.

A cooler body improves the ability to sustain intensive workloads over longer periods of time. The THORZT Chill towel is a great way to cool the skin and contributes to a number of physiological advantages:

- Delay the onset of dehydration
- Less output is directed toward the skin
- A lower core body temperature

The THORZT Chill Towel stays cool for hours with no refrigeration or cold water needed - just 'Wet, Wring and Wave'. It's made from PVA material which can retain moisture for periods far longer than conventional towels, and contains a special plant extract that (once activated) turns moisture cold regardless of it's initial temperature.

FEATURES & BENEFITS

- Promotes Thermoregulation
- Reduces Heat Stress risks
- Comforts and cools the skin
- Looks as cool as it feels
- Machine washable

PRODUCT DETAILS

Size: 75cm x 32cm

Ctn Qty: 20

CARE INSTRUCTIONS

- It is recommended that the Chill Towel is washed or rinsed before use.
- When finished using, fold the damp Chill Towel and return to packaging - do not fold when dry.
- Machine wash with mild detergent.
- Do not put in dryer.



SCAN TO
WATCH VIDEO



CSB

CSP

ACTIVATION

1. Soak and Submerge the THORZT Chill Towel in any temperature water.
2. Wring out the excess water, leaving it nice and moist.
3. Take hold of two corners, and wave it up and down for about 10 seconds.
4. Wrap it around your neck or other warm point on the skin for instant cooling relief.

REACTIVATION

Give it a quick wave every 15-20 minutes to reactivate the cooling process. Re-wet every few hours if needed.